
## ZAPAMTITE

## Metode za ublažavanje interferencije daratumumaba

Kod pacijenata lečenih daratumumabom moguća je paneaktivnost kod indirektnog antiglobulinskog testa (IAT)

Metode za ublažavanje interferencije daratumumaba

III

Ako je dostupna, pogledati identifikacionu karticu za pacijente koja sadrži podatke o njihovoj krvnoj grupi i rezultata skrininga na antitela urađenog pre početka lečenja daratumumabom.

## DARATUMUMAB

Važne informacije  
o minimizaciji rizika od interferencije  
daratumumaba sa testovima  
kompatibilnosti krvi!  
Brošura za hematološke / transfuzijske  
laboratorije

## REFERENCE:

1. Chapuy CI, Nicholson RT, Aguad MD, et al. Resolving the daratumumab interference with blood compatibility testing. Transfusion. 2015;55(6 Pt 2):1545-155.
2. Albeniz I, Demir O, Türker-Sener L, Yalcintepe L, Nurten R, Bermek E. Erythrocyte CD38 as a prognostic marker in cancer. Hematology. 2007;12(5):409-414.
3. Mehta K, Shahid U, Malavasi F. Human CD38, a cell-surface protein with multiple functions. FASEB J. 1996;10(12):1408-1417.
4. Zocchi E, Franco L, Guida L, et al. A single protein immunologically identified as CD38 displays NAD+ glycohydrolase, ADP-ribosyl cyclase and cyclic ADP-ribose hydrolase activities at the outer surface of human erythrocytes. Biochem Biophys Res Commun. 1993;196(3):1459-1465.
5. Oostendorp M, Lammerts van Bueren JJ, Doshi P, et al. When blood transfusion medicine becomes complicated due to interference by monoclonal antibody therapy. Transfusion. 2015;55(6 Pt 2):1555-1562.
6. Hannon JL, Clarke G. Transfusion management of patients receiving daratumumab therapy for advanced plasma cell myeloma. Transfusion. 2015;55(11):2770.
7. Westhoff CM, Reid ME. Review: the Kell, Duffy, and Kidd blood group systems. Immunohematology. 2004;20(1):37-49.
8. Sažetak karakteristika leka Darzalex, dostupan na internet stranicama <https://cinmed.me/>

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